

# Microstructural Analysis of Solar Cells

[Download PDF Copy](#)[View Supplier Profile](#)[Request Quote](#)Sponsored by [Buehler](#)

Apr 14 2016

## Table of Contents

[Introduction](#)[Common Anomalies](#)[Sampling Techniques](#)[Sample Mounting](#)[Sample Grinding](#)[Sample Polishing](#)[Final Polishing](#)[Metallographic Examination](#)[Silicon Ingot Sample Preparation](#)

## Introduction

The development and commercialization of photovoltaic (PV) cells is being carried out rapidly. Manufacturers of wafer cells and modules are expected to manufacture panels that are extremely reliable and durable. Although there is a huge demand for warranties of more than 25 years in the market, the quality of a product can deteriorate over time due to the persistent stress on the metallic interconnects. This issue can be overcome by performing [cross-sectional microscopic analysis](#) to enhance module/cell reliability from material inspection and interconnect verification of coating analysis. This article presents different techniques to prepare high-quality cross-section samples of crystalline silicon solar cells to

This site uses cookies. By continuing to browse the site you are agreeing to our use of cookies. [Find out more.](#)




## Common Anomalies

Exposure to thermal cycling is the main cause of failures in electronic assemblies. To resist repeated thermal stress, tab ribbon bonds should be durable and reliable. To assess thermo-mechanical stress during pre- and post-accelerated aging of cell interconnects, microstructural analysis is carried out. Aging of the solder interconnects causes constant expansion and contraction of the joint, which causes a solder to fatigue, become fragile, and break into larger grains, altering the mechanical properties of the solder interconnects. To evaluate the efficiency of old and new bonding materials, or when transitioning from leaded to lead-free solders, cell research is generally carried out. Cross sectional analysis of raw materials, microelectronic electric components, and wafers reveal the parameters mentioned below:

- Poor solder wetting
- Joint cracks
- Impurities
- Solder porosity and microvoids
- Solder thickness and delamination
- Solder meniscus
- Unacceptable solder, silicon, and copper microstructure
- Intermetallic phases, coarsening, etc.

## Sampling Techniques

A solder joint, copper ribbon, and bus bar (metalized lines) can be extracted from a cell to prepare micro-sections. Solar cells are serially connected on back and front sides (Figure 1) with copper solder coated ribbons. These ribbons are soft soldered (lead based solder) to the paste-like bus bar material, which mainly consists of glass and silver particles. The latter enables solderability and electrical conductivity for the cells. To bond the bus bar material to the cell before the application of the solder-coated ribbons (Figure 2), it is hardened in an oven. To prolong the service life and enhance the electrical conductivity, the use of high-quality solder joints are essential.




**Figure 1.** Solar cells connected in series with tabbing ribbons



**Figure 2.** Interconnector on a soldered solar cell sample exhibiting a poor quality cut (top) and a good quality cut (bottom)

Sectioning is carried out on monocrystalline (thick-films) or polycrystalline solar cells to downsize the sample, and bring it to a manageable size to mount, grind and polish. To carry out cross-sectional analysis of a PV cell interconnect, the wafer and the solder joint (Figure 3) have to be sectioned without causing any major damage to the cut surface. This preparation is tedious, as the materials are fragile and possess a wide variety of mechanical properties. It is recommended to perform encapsulation with a quick drying acrylic resin such as SamplKwick®, which cures in 15 to 25 minutes to inhibit chipping and cracking during sectioning. The resin is applied by painting the bus lines onto the sides of the cell with a brush (Figure 4). The resin may need to be removed after cutting and before potting to render a stronger bond of the sample while being mounted in epoxy.




**Figure 3.** Solder interconnect of a solar cell bus line



**Figure 4.** *Front side, top and back side of a polycrystalline Silicon solar cell*

After the acrylic resin covering the bus lines is completely cured, the whole cell, measuring 156 x 156 mm x roughly 0.2 mm in size, is positioned in a solar cell clamping device ([Buehler #11-2706 holder](#)) employed in the [IsoMet<sup>®</sup> 4000 Linear Precision Saw](#) (Figure 5). It should be noted that when utilizing a solar cell holder, an acrylic resin coating may no longer be needed.



**Figure 5.** Buehler solar cell holder for IsoMet® 4000 and 5000 Linear Precision Saws


After sectioning the entire bus bar, it is cross cut to a manageable size for subsequent mounting (Figure 6). Samples measuring roughly 1 cm<sup>2</sup> are cut from the cell at a cutting speed of 3750 rpm and a feed rate of 5 mm/minute using an IsoMet® Diamond Wafering Blade LC-15 type. The suitable selection of blades reduces the damage that would have to be eliminated by further grinding and polishing. Typical issues experienced while carrying out sectioning are; silicon breakage, bus bar detachment, and edge cracking. Cracks or broken edges (Figure 2) occur due to an extremely high feed rate or a coarser blade.



**Figure 6.** Sectioning of the bus bar into samples appropriate for mounting

## Sample Mounting

After the sample has been sectioned from the PV solar cell, it is mounted in epoxy by casting in an epoxy resin EpoxiCure® (Figure 7). EpoxiCure® ensures optimum sample edge retention and is cured within 6 hours, without any notable increase in the temperature in the mounting cup measuring 1.25" (32 mm). EpoxiCure® can be used in combination with a Conductive Filler to perform energy dispersive X-ray spectrometer (EDS) analysis, or scanning electron microscopy (SEM). For vertical positioning of the sample in the SamplKup® mounting cup, a UniClip Support Clip (Figure 8) is utilized. In case of additional edge support requirements, the epoxy is combined with a Flat Edge Filler (Figure 9). The Flat Edge Filler consists of globular ceramic oxide particles, which increase the mount hardness and minimize shrinkage. Sample edge retention is a key feature when observing features at the edge, or when measuring the coating thickness.



**Figure 7.** Solar cell interconnects mounted in EpoxiCure® epoxy



**Figure 8.** Solar cell with interconnector held perpendicular to the cross section with a UniClip Support Clip before epoxy mounting.



**Figure 9.** EpoxiCure<sup>®</sup> mount solar cells interconnect with Flat Edge Filler (white layer) added to improve edge retention.

## Sample Grinding

After being mounted the samples are mechanically ground and polished. Initially, planar grinding of the samples is carried out using the CarbiMet<sup>®</sup> 2, 320 grit sized SiC grinding paper. Grinding eliminates any sectioning damage and ensures a flat surface. As induced artifacts can appear upon grinding, it is recommended to maintain a minimal grinding time. The ductile and soft metals (like solder and copper) that trap the abrasives break off upon grinding, and are usually removed at the time of polishing. With each preparation step after grinding care should be taken to thoroughly clean the sample holder that is utilized with semi-automated grinder polishers, and to clean machine components and samples with tap water. When using an UltraMet<sup>®</sup> Ultrasonic Cleaner, distilled water should be used to clean the samples, in order to prevent scratches in the sample when some of the abrasives are carried to the next step. Scratches can cause problems with electron backscatter diffraction (EBSD) microstructural analysis and fault analysis.


## Sample Polishing

Following grinding, the samples are refined by multiple polishing steps. Polishing of Si cells involves four steps, including polishing with alumina and diamond abrasives charged on a polishing cloth. A finer abrasive, lubricated with MetaDi<sup>®</sup> Fluid, is utilized in each step. As a result of varied hardness of the cell material, hard polishing cloths are used to obtain excellent results. When compared with silicon which is hard, copper and solder are relatively soft and ductile. When the hard and soft materials are placed near each other in a solar cell the soft materials are removed at a greater pace, resulting in excessive relief. This makes sample analysis complicated, and reduces the depth of field because of greater differences in material heights.


## Final Polishing

Final polish is the last step of the standard technique, where extremely fine abrasive (0.02  $\mu\text{m}$ ) silica are used to eliminate fine scratches remaining from the prior step. Only a small amount of material is removed, so any damage that has been overlooked might not be rectified. Final polishing should be completed within 2 minutes if there has been proper preparation of the sample before this step. The sample can undergo a final preparation vibratory polishing at this point to eliminate any smearing effect, and to clearly define between materials and certain alloy phases. For copper and solder alloys, scratch removal may be extremely hard. To eliminate slight scratches and surplus surface damage, vibratory polishing is performed with a low nap polishing cloth and a colloidal silica suspension. Figures 10 and 11

illustrate the outcome of using a vibratory polisher on an SAC 305 solder alloy. The entire sample preparation method and related parameters are provided in Table 1.



*Figure 10. Unetched SAC 305 solder alloy*



*Figure 11. Unetched SAC 305 Solder Photomicrograph after vibratory polishing*


*Table 1. Sample preparation of polycrystalline solar cell solder interconnect*




Step	Surface	Lubricant	Abrasive/ Size	Time (mins.)	Force Lb (N)	Base/ Head rpm	Head Direction
1. Surface Grind	CarbiMet <sup>®</sup> 2 SiC Paper	Water	320 (P400) grit	Until Plane	5 lb (22)	300/60	Complementary
2.	UltraPad <sup>™</sup>	MetaDi <sup>®</sup>	9 μm MetaDi <sup>®</sup> Supreme Diamond	5:00	5 lb (22)	150/60	Contra
3.	TexMet <sup>™</sup>	MetaDi <sup>®</sup>	3 μm MetaDi <sup>®</sup> Supreme Diamond	4:00	5 lb (22)	150/60	Complementary
4.	MicroCloth <sup>®</sup>	MetaDi <sup>®</sup>	1 μm MetaDi <sup>®</sup> Supreme Diamond	2:00	5 lb (22)	150/60	Complementary
5.	ChemoMet <sup>®</sup>	MetaDi <sup>®</sup>	0.02 μm MasterMet <sup>®</sup> 2 Colloidal Silica Suspension	1:30	5 lb (22)	150/60	Contra
6. Optional	VibroMet <sup>®</sup> 2 Vibratory Polisher; MicroCloth <sup>®</sup>	N/A	0.02 μm MasterMet <sup>®</sup> 2 Colloidal Silica Suspension	60	N/A	N/A	N/A

## Metallographic Examination

The quality of solar modular production is ensured by investigating the soft solder joints. Solder interconnect quality plays a vital role in increasing the lifetime of solar cells. Individual cells are serially connected by soldering flat copper solder coated ribbon to metalized bar traces. Figure 12 depicts a cross-sectional view of the ribbon interconnect. Figure 13 shows insufficient solder at each end of the solder encapsulated copper ribbon. Measurement can be made to ensure proper thickness of the other solder layers and bus bar.




**Figure 12.** Microstructure of Sn62-Pb36-Ag02 solder, unetched solar cell copper ribbon




**Figure 13.** The entire solar cell interconnect shows both top and bottom soldered seams.

Wetting is an alloy process where the solder coated ribbon is soldered to the bus bar. To ensure strong stable connections it is crucial to perform proper wetting. In the field of electronics it is possible to determine a wetting angle, and a slightly concave wetting angle of less than  $40^\circ$  is sufficient. As screen-printed pastes are used for the bus lines in solar cells it is hard to attain this angle. A satisfactory solder wetting angle between the soldered coated copper ribbon and the bus bar is shown in Figure 14. Figure 15 displays an inappropriate wetting angle showing voids in the solder. Depending on their size and location, a specific percentage (area) of voids is allowed. Microvoiding (small voids along the interface), usually less than 1 mm in diameter, are often caused due to the surface finish or contamination.



**Figure 14.** Proper wetting angle of soldered copper



**Figure 15.** Both solder interconnects show an insufficient wetting angle and several voids in the solder layer

Thermal stresses can lead to delamination and the formation of cracks during wafer/cell manufacturing, or during soldering. So when preparing the sample, care must be taken to avoid smearing of soft metals that hide micro-cracks. Figure 16 displays a delamination of the bus bar and a crack in the silicon wafer.



**Figure 16.** Light optical images showing delamination of a contacted bus bar at the cell interconnect and cracking seen across the wafer in dark gray at the bottom


When two dissimilar metals diffuse into one another, intermetallic compounds can form. Numerous lead-free and leaded solders will develop over the lifetime of the cell. These compounds may become hard and brittle, and they can have greater electrical resistivity than the constituent metal. A bismuth with a lead-free solder interconnect is shown in Figure 17. Figure 18 shows Sn62-Pb36-Ag02, etched with Klemm II, Eutectic solder material, with Pb-rich solid solution and Sn-rich solid solution. Figure 19 displays SEM images of the intermetallic phase growth of the solder above the bus bar, both before and after the accelerated aging of a solder interconnect on a bismuth lead-free solder. The different compounds developed in this process are analyzed with EDX by labeling the materials and their elemental maps. Interpretation of the constantly evolving microstructure and its mechanical properties when subjected to accelerated aging, including the effects on the solder joint, would provide a better understanding about how to increase the lifetime of the interconnects.



**Figure 17.** Microstructure of Bismuth containing lead-free solder Sn60-Bi38-Ag02 (unetched) interconnect



**Figure 18.** Sn62-Pb36-Ag02, etched with Klemm II, Eutectic solder material with Pb-rich solid solution (dark) and Sn-rich solid solution (light-colored)



**Figure 19.** SEM image showing the intermetallic phase growth after accelerated aging: (a) Bismuth containing solder on (b) bus bar

## Silicon Ingot Sample Preparation

Performing microstructural analysis of amorphous silicon ingots, utilized in wafer cell production, is crucial to enhance product performance. The characterization of the material for quality purposes can expose grain size, structure, impurities, annealing twins, and give a record of its process for future purposes. Preparation of pure silicon or pure metals is a complicated and time-consuming process compared to alloys.

This process starts with removal of the material from the ingot, but most of the damage occurs during sectioning, so care must be taken to ensure that the right sectioning blade is used to reduce thermal and mechanical damage. It is advisable to perform compression mounting using a [SimpliMet<sup>®</sup> automated Mounting Press](#) when the sample is ready to be mounted, as it is not as delicate as the solar cell interconnects. To ensure better hardness and edge support, [EpoMet<sup>®</sup> F Molding Compound](#) is utilized for the mount media. After mounting the sample, the sample preparation can be performed with a seven-step preparation method (Table 2).


**Table 2.** Pure silicon preparation method

Step	Surface	Lubricant	Abrasive/ Size	Time (mins.)	Force lb [N]	Base/ Head rpm	Head Direction
1. Surface Grind	CarbiMet <sup>®</sup> 2 SiC Paper	Water Cooled	240 grit [P280] SiC	Until Plane	5lb [22]	300	Complementary
2.	CarbiMet <sup>®</sup> 2	Water	320 grit [P400] SiC	1	5lb [22]	300	Complementary


	SiC Paper	Cooled					
3.	UltraPad™ Silk	MetaDi® Fluid	9 μm MetaDi®	5	5lb [22]	150	Contra
4.	TriDent®	MetaDi® Fluid	3 μm MetaDi®	5	5lb [22]	150	Complementary
5.	TriDent®	MetaDi® Fluid	1 μm MetaDi®	2	5lb [22]	150	Complementary
6.	MicroCloth®		0.02 μm MasterMet® 2 Colloidal Silica	3	5lb [22]	150	Contra
7.	MicroCloth®		0.02 μm MasterMet® 2 Colloidal Silica	60	VibroMet® 2 Vibratory Polisher		




**Figure 20.** High purity, defect-free, polycrystalline silicon viewed with Nomarski DIC (etched with 100mL water, 75g NaOH)



**Figure 21.** Polycrystalline silicon viewed with Nomarski DIC (etched with 100mL water, 75g NaOH)



**Figure 22.** Standard sample preparation of pure single-crystal, band contrast 205.8



**Figure 23.** Standard sample preparation plus vibratory polishing of pure single-crystal, band contrast 233

## References

1. Schmitt P, Eberlein D, Voos P, Tranitz M, Wirth H. Metallographic Preparation of Solar Cell Samples for Quality Assurance and Material Evaluation. 1st International Conference on Silicon Photovoltaics, April 17-20, 2011, Germany, 2011.
2. Practical Metallography: Eberlein D, Schmitt P, Voos P. Metallographic Sample Preparation of Soldered Solar Cells. Practical Metallography 2011; pp. 239 – 260; 5-2011; 2011.
3. G. Cuddalopepatta, et al.: Durability of Pb free solder between copper interconnect and silicon in photovoltaics cells. in: Prog. Photovolt. Res. Appl., Vol. 18, pp. 168-182.
4. BUEHLER SUM-MET – The Science Behind Materials Preparation, BUEHLER LTD, 2007.
5. George F. Vander Voort: ASM Handbook Volume 9 – Metallography and Microstructures, 2004.



This information has been sourced, reviewed and adapted from materials provided by Buehler.

For more information on this source, please visit [Buehler](#).

Download PDF Copy

View Supplier Profile

Request Quote

Read in: English

## Ask A Question

Do you have a question you'd like to ask regarding this article?



Leave your feedback

Login    ... 

Public Question  Private Question to Supplier

Submit

Helping engineers find the right adhesive system

- ✓ EPOXIES
- ✓ SILICONES
- ✓ UV CURES



 **MASTERBOND**

[www.masterbond.com](http://www.masterbond.com)



Latest interviews

Featured equipment

Trending stories

### Cooling Generators with Hydrogen

In this interview, Thomas Skoczylas, Territory Sales Manager from Proton OnSite talks to AZoM about cooling generators with hydrogen.




### Optical Metrology and its Current Trends

Optical Metrology has found increased demands for difficult environments and increased data storing capabilities. AZoM spoke Dr. Peter de Groot, of ZYGO, about his thoughts and expertise surrounding the subject and how he feels ZYGO will lead the industry moving forward.



### Accurate and Aesthetic Viscosity – Introducing the V-Compact

Fungilab’s main focus since its establishment has been to research, develop and manufacture the most advanced viscometers in order to offer our clients the highest and most accurate performance in the viscosity measurement. In this interview, AZoM speaks to Fungilab CEO, Ernest Buira.



**OLYMPUS**




DSX Digital Microscopes

**OLYMPUS**

Your Turnkey Solution  
for Technical  
Cleanliness Inspection

**The CIX100 System**



**KEYSIGHT TECHNOLOGIES**



**SHIMADZU**

High-performance  
UV-Vis-NIR

Routine, accurate  
characterization  
of a variety of  
materials

## Related ANSI Standards

**ASTM B193-16:** Standard Test Method for Resistivity of Electrical Conductor Materials

**ASTM B593-96(2014)e1:** Standard Test Method for Bending Fatigue Testing for Copper-Alloy Spring Materials

**ASTM B702-93(2015):** Standard Specification for Copper-Tungsten Electrical Contact Material


**ASTM B795-13:** Standard Test Method for Determining the Percentage of Alloyed or Unalloyed Iron Contamination Present in Powder Forged (PF) Steel Materials

**ASTM E110-14:** Standard Test Method for Rockwell and Brinell Hardness of Metallic Materials by Portable Hardness Testers


## More Content from Buehler



Maximizing the Potential of Your Abrasive Saw



Fastener Metallography: Assuring the Highest Quality of Compliance



Resin Bonded Diamond Surface Enhance Microstructural Analysis of Thermally Sprayed Coatings

See all content from Buehler



Find, evaluate and source the best materials for your projects.

## Other Sites from AZoNetwork

AZoNano  
AZoRobotics  
AZoCleantech  
AZoOptics  
AZoBuild

AZoMining  
AZoSensors  
AZoQuantum  
News Medical  
Life Sciences

## Useful Links

News  
Articles  
Directory  
Equipment  
Books  
Videos  
Events  
Metals Store  
Journals  
Software

Courses  
Webinars  
Events  
Market Reports  
Materials  
Applications  
Industries  
AZojomo  
Podcasts  
Experts  
The Team

Search  
Become a Member  
Newsletters  
About  
Contact  
Help/FAQs  
Advertise  
Terms  
Sitemap

[back to top](#)



AZoM.com - An AZoNetwork Site

Owned and operated by AZoNetwork, © 2000-2017