






BUEHLER® Sum-Met™ Bundles

- 1 Part Number
- 1 Price
- Follows Sum-Met™ Preparation Guide
- Each Bundle Prepares 100 Specimen Holder
- Preprogrammed in EcoMet® Pro Grinder/Polisher



Material	Bundle Contents	8" (203mm)	10" (254mm)	12" (305mm)
Fe & Ni	 CarbiMet® 2 (320/P400), UltraPad™, Trident™ & MicroCloth® polishing Cloths, 9µm and 3µm MetaDi® Supreme Diamond Suspension, MasterPrep® Final Polishing Suspension	42-10017	42-10018	42-10019
Al & Cu	 CarbiMet® 2 (320/P400), UltraPad™, TexMet® & MasterTex® polishing Cloths, 9µm, 3µm and 1µm MetaDi® Supreme Diamond Paste, MetaDi® Fluid and MasterMet® Final polishing Suspension	42-10020	42-10021	42-10022
Ti	 CarbiMet® 2 (240/P280), UltraPol™ and ChemoMet® Polishing Cloths, 9µm MetaDi® Supreme Diamond Suspension, MasterMet® Final Polishing Suspension	42-10023	42-10024	42-10025
Electronic & Polymer	 CarbiMet® 2 (240/P280) and (600/P1200), TexMet® and MicroCloth® Polishing Cloths, 3µm MetaDi® Supreme Diamond Paste, MetaDi® Fluid and MasterPrep® Final Polishing Suspension	42-10026	42-10027	42-10028
Ceramic & Refractory	 Apex® DGD Red Disc, UltraPol™, TexMet® and ChemoMet® Polishing Cloths, 9µm and 3µm MetaDi® Diamond Supreme and MasterMet Final Polishing Suspension	42-10029	42-10030	42-10031

For a complete listing of consumables, please refer to our Consumables Buyer's Guide.
Buehler continuously makes product improvements; therefore, technical specifications are subject to change without notice.

BUEHLER GMBH - European and MESA Headquarters

In der Steele 2 • 40599 Düsseldorf
Postfach 16 03 55 • 40566 Düsseldorf
Telefon: (49) 211 974100 • Telefax: (49) 211 9741079
www.buehler-met.de
info@buehler-met.de

BUEHLER FRANCE
Téléphone: 0800 89 73 71
Télécopie: 0800 88 05 27
www.buehler.fr
info@buehler.fr

BUEHLER UNITED KINGDOM
Telephone: 0800 707 6273
Fax: 0800 707 6274
www.buehler.co.uk
sales@buehler.co.uk

BUEHLER CANADA

10 Carlow Court, Unit #2
Whitby, Ontario L1N 9T7
Telephone: (905) 430-4684 • Fax: (905) 430-4647
Sales Telephone: 1-800-268-3593
info@buehler.ca

BUEHLER, ASIA-PACIFIC

5/F Vogue Centre
696 Castle Peak Road
Lai Chi Kok, Kowloon
Hong Kong, SAR, China
Telephone: (852) 2307 0909
Fax: (852) 2307 0233