



BUEHLER

Buehler, Ltd.

41 Waukegan Road
Lake Bluff, IL 60044
www.buehler.com
Phone: 847-295-6500
Fax: 847-295-7979
E-mail: info@buehler.com

PRESS RELEASE

For Immediate Release

Contact: Peter Leal

E-mail: Peter.Leal@buehler.com

Buehler's New AbrasiMatic® 300 Abrasive Cutters!



Lake Bluff, Il., July 15th, 2009- Buehler, Ltd. proudly introduces our new AbrasiMatic® 300 family of abrasive cutters!

Buehler's new AbrasiMatic® 300 Abrasive Cutter is a versatile bench-top cutter featuring manual cutting action in three directions and automated cutting action in one direction. Features include a 12" (305mm) cut-off wheel capacity, high torque motor and large cutting capacity. Buehler's unique SMARTCUT system automatically reduces feed rate to produce burn-free, consistent sections. The user friendly color touch-screen control provides the maximum versatility to section a wide variety of sample materials, sizes and geometries!

#

Buehler Ltd. is the world's premier manufacturer of scientific equipment and supplies for use in materials analysis. Buehler's products are used throughout the world in manufacturing facilities, quality laboratories, and universities to enable material characterization, ensure quality, and perform materials research. For more information, contact Buehler, Ltd., 41 Waukegan Road, Lake Bluff, IL 60044, (ph) 847-295-6500.