



BUEHLER

Buehler, Ltd.

41 Waukegan Road
Lake Bluff, IL 60044
www.buehler.com
Phone: 847-295-6500
Fax: 847-295-7979
E-mail: info@buehler.com

PRESS RELEASE

For Immediate Release

Contact: Peter Leal

E-mail: Peter.Leal@buehler.com

The Future of Material Sample Preparation is Here!



Lake Bluff, Il., April 16th, 2009- Buehler, Ltd. proudly introduces the new EcoMet® 300 Grinder-Polisher and AutoMet® 300 Power Head, Buehler's newest family of robust, high performing metallographic grinder-polishers.

The EcoMet® 300 Grinder-Polisher and AutoMet® 300 Power Head are designed for the materials analyst who requires versatility in both manual and semi-automatic sample preparation with consistent and repeatable results. The EcoMet® 300 Grinder-Polisher can be used alone for manual grinding applications or combined with the AutoMet® 300 Power Head for semi-automatic applications. Features include variable speed controls, membrane keypad, unique LED lighting, retractable water hose, quick release chuck, replaceable bowl liner, and a 360° bowl rinse system. Users can now successfully prepare anywhere from a single simple sample to multiple, unique geometry shaped samples.

###

Buehler is the world's premier manufacturer of scientific equipment and supplies for use in materials analysis. Buehler's products are used throughout the world in manufacturing facilities, quality laboratories, and universities to enable material characterization, ensure quality, and perform materials research. For more information, contact Buehler, Ltd., 41 Waukegan Road, Lake Bluff, IL 60044, (ph) 847-295-6500.