




# Premium Resin Abrasive Blades

*Longer Life,  
Great Quality*

Buehler's new resin abrasive blades are individually targeted to cut a full range of metallic materials and last up to 3 times longer than competitor blades. Users can rely on the consistency of the blades while the longevity leads to fewer blade changeovers and purchases.




## Efficient

Buehler's premium resin abrasive blades cut through samples fast without sacrificing cut quality. Have confidence in the consistency, cut quality and new prolonged life of these blades to improve your preparation process.



## QR Scanning

Buehler's new abrasive blades come with QR codes to access blade information and simplify the reordering process. Simply scan and follow the link to place your order.




This line of resin blades is compatible with all current Buehler equipment and any machine with a 1.25in [32mm] arbor.

The resin construction ensures users will not be bothered by the excessive odor that rubber bonded blades can produce.

## Ordering Information

### 1.25in [32mm] Arbor (Qty 10)

Recommended Use	10in [254mm]	12in [305mm]	16in [406mm]
Non Ferrous Materials			10-31650-010
Ferrous Materials HRC 35-50	10-31014-010	10-31218-010	10-31616-010
Ferrous Materials HRC 50-60	10-31014-010	10-31212-010	10-31612-010
Ferrous Materials >HRC 60		10-31218-010	10-31610-010
Super Alloys		10-31205-010	10-31605-010



### All New

The 16in [406mm] abrasive blades are perfect to use with Buehler's new AbrasiMet XL Pro (Item 10-2320-400). This floor standing cutter, which can accommodate blades up to 18in [457mm], meets all of your abrasive cutting needs. The 10in [254mm] and 12 in [305mm] abrasive blades are perfect to use with Buehler's new AbrasiMet M (Item 10-10107-200/400). This benchtop cutter is great for cutting small and medium sized samples at a great value. Contact your sales representative to request a quote! <https://bhlr.io/2Z4MqhG>



**BUEHLER Worldwide Headquarters  
North America-South America Offices**  
41 Waukegan Road  
Lake Bluff, Illinois 60044-1699 USA  
P: 800 BUEHLER (800-283-4537)  
P: (847) 295-6500  
[www.buehler.com](http://www.buehler.com) | [info@buehler.com](mailto:info@buehler.com)

**European Headquarters  
BUEHLER Germany**  
[info.eu@buehler.com](mailto:info.eu@buehler.com)

**BUEHLER France**  
[info.eu@buehler.com](mailto:info.eu@buehler.com)

**BUEHLER United Kingdom**  
[info.eu@buehler.com](mailto:info.eu@buehler.com)

**BUEHLER China**  
[info.cn@buehler.com](mailto:info.cn@buehler.com)

**BUEHLER Japan**  
[info.japan@buehler.com](mailto:info.japan@buehler.com)

**BUEHLER Asia-Pacific**  
[info.asia@buehler.com](mailto:info.asia@buehler.com)