



PREPARE LARGE, COMPLEX SAMPLES WITH CONFIDENCE

Preparing large or unmounted specimens for semi-automatic grinding & polishing is time consuming and difficult. Buehler's Universal Sample Holder solve these challenges and simplify your process.

Secure Samples, High Quality Results

Universal Universal Sample Holders securely clamp large and irregular specimens using adjustable supports and clamping pins. This ensures stable positioning, consistent force, and improved safety and preparation quality.

Flexible Configuration for Your Lab

The modular design adapts to a wide range of sample sizes, shapes, and materials. Adjustable support columns and optional polymer-coated pins allow secure clamping, even for brittle samples.



Ordering Information

Available in 127 mm, 157 mm, and 200 mm widths to match different platens and sample sizes. Each holder includes support columns and metal-tipped pins for a customizable setup. Polymer-tipped screws are available separately.

Part Number	Description
492505	127mm Universal Sample Holder
492506	157mm Universal Sample Holder
492507	200mm Universal Sample Holder
492508	Polymer Coated Pins

