



Scan to  
Learn about  
Abrasive  
Vising

## FAST. ROBUST. EFFORTLESS.

Boost your workflow with this efficient vise designed to make clamping fast and easy. Its streamlined build and drop-in tightening screw allow quick adjustments, helping you save time and work with precision.


### Technical Specifications

**Size:** Medium

**Max Opening:** 2.75in [69.8mm]

**Clamping Plate:** 2in x 1.5in [50mm x 40mm]

**T-Nuts:** 12mm




**Part Number: 10-3555**

### Compatible with:



**IsoMet®  
High Speed  
Precision Cutter**



**AbrasiMet® M  
Abrasive Cutter**



**AbrasiMatic® 300  
Abrasive Cutter**

