

# WIRE AND TUBE NEWS

**MAC**  
Magnetic Analysis Corp.

**Multi-Test System  
For Shock Absorber Tubing**

HOME WIRE TUBE WIRE EVENTS TUBE EVENTS CONTACT US

SEARCH

Events - Tube Events - Wire

## Buehler at Control 2016: From Sample Preparation to Automated Hardness Testing

📅 July 2, 2016

At Control 2016, Buehler ITW Test and Measurement will be showcasing their latest additions to their range of products for metallography preparation and hardness testing, designed for applications ranging from fundamental research to automated, in-process quality control testing.



Buehler ITW Test & Measurement ([www.buehler.com](http://www.buehler.com)) will be introducing, at its booth no. 5135 in Hall 5, new solutions from its extensive product portfolio, comprising equipment, consumables and accessories for metallography preparation and

hardness testing, all under one roof. For the first time, Buehler will be showcasing its ergonomically optimized VH1002 Series hardness testers – configurable in various levels of automation and equipped with weights from 10 g up to 2 kg – as well as its newly developed Diamet™ universal hardness testing software.

Another highlight at the booth will be the new SimpliMet™ 4000 mounting press – a highly reliable and time-saving solution for mounting of samples prior to grinding and polishing, specifically designed for use in a demanding industrial environment.

Further exhibits will include the recently launched PlanarMet™ 300 table-top grinder and the Wilson® VH3300 micro hardness tester, both of which provide for a high degree of automation. Exhibits from the areas of grinding, polishing, sectioning



and mounting of samples illustrate the width of Buehler's portfolio, which covers applications from fundamental research to automated in-process quality control.

The revolutionary DiaMet™ software suite is designed for analyzing macro and micro Vickers, Knoop, Rockwell and Brinell indents in accordance with all relevant ISO and ASTM standards. The software is the first solution which can be used across a wide range of systems and is tailored to meet global hardness testing requirements, enabling measurements to be completed and documented reliably and compared internationally across a wide range of industries. Featuring clear menu navigation and high ease of operation, DiaMet™ meets the requirements of a variety of users, from production staff with limited training up to specialists in scientific research.

PlanarMet™ 300 is a recently released bench-top planar grinder which provides rapid, highly accurate and at the same time cost-efficient and convenient sample preparation in metallography. The unit's fixed 305 mm grinding stone along with a high-performance motor is designed to remove up to 0.6 mm of material per minute – up to ten times faster than grinding with SiC paper. The central-force specimen holder, accommodating up to 10 samples at a time, enables three to four grinding stages to be combined into a single step. The time and cost saving PlanarMet™ 300 system therefore offers all the power of a large, floor-standing model in a space-saving, compact bench-top design.

The Wilson® VH3300 Vickers and Knoop hardness testers were designed for maximum flexibility and efficiency. These systems excel through a wide range of test loads from 10 g to 50 kg, a high magnification and a 6-position test head with three indenters and three zoom objectives, all of which contribute to its universal applicability. A patent-pending collision protection system and a high-speed sample stage system contribute towards speeding up the workflow in automated hardness testing. The system's area of application extends from material and product development to routine quality control testing in a clean production environment.

Buehler ITW Test & Measurement's product portfolio includes a wide variety of sectioning and high-precision sectioning systems, mounting systems including the associated epoxy and acrylic resins, grinders and polishers as well as Rockwell, Vickers/Knoop, Brinell and universal hardness testers.

**Buehler – ITW Test & Measurement GmbH**, based in Esslingen, Germany has been a leading manufacturer of scientific equipment, supplies and accessories for use in metallography and materials analysis for 80 years now, and also offers a broad range of hardness testers and hardness testing systems. A close-knit network of branch offices and dealers provides customers with expert support and service around the globe. The Buehler Solutions Centre in Esslingen/Germany as well as many other Solutions Centres across Europe and the world offer extensive assistance with all application-related questions and the development of repeatable specimen preparation processes.

## TUBE NEWS

October 2017

60 Years  
of ZSK Twin  
Screw Extruders  
Anniversary of Coperton's  
high-performance ZSK  
extruder series



### INSIDE

Tube Processing on a  
Flat-Bed Fibre Laser  
Cutting Machine

Highest Quality Welds  
made with Lower Costs

Ring Roller Closes the  
Gap

## Magazine Archive

🔖 TN 2015

🔖 TN 2016

🔖 TN 2017

🔖 WN 2015

🔖 WN 2016

🔖 WN 2017

**Buehler** is part of the Test and Measurement segment of the US based Illinois Tool Works (ITW) group of companies with more than 200 distributed business units in 52 countries worldwide, and a staff of approx. 45,000.

Contact: +49 (0) 711 4904690-847

Share this:



#### Related



[Buehler at Control 2017](#)

In "Events - Wire"



[EpoKwick™ FC Fast Cure Epoxy Mounting Compound for Top Quality Sample Preparation in Materialography](#)

In "Events - Tube"



[New SimpliMet™ 4000 Mounting Press Shortens Cycle Times and Increases Reliability in the Preparation of Metallographic Samples](#)

In "Wire"

[← AICON 3D Systems](#)

[Metalube Expands Manufacturing Capacity →](#)