



Delta Manual Service Plan



Protect Your Investment

Labs run more efficiently and consistently with well-maintained equipment.



Reduce Unplanned Downtime

Regular equipment maintenance will reduce the hassle of unexpected downtime and associated costs.



Ensure Consistent Results

Confidently analyze lab results with the equipment maintained by the trusted experts at Buehler.



Preventative Maintenance

Protect your investment and ensure consistent performance, minimize downtime and reduce the costly repairs with Buehler Preventative Maintenance and Calibration Services.

Preventative Maintenance Service Plan Ensures Consistent and Reliable Results

- Annual visit includes Wear Parts Kit and 15 Point Maintenance Check list
- Trained technicians able to perform Preventative Maintenance on multiple machines in one visit
- Unlimited Phone Support
- 10% Replacement Part Discount





Maintenance Checklist

Our comprehensive 15 point maintenance checklist is thorough, ensuring the Delta Manual is running efficiently.

15 Point Maintenance Checklist

- Check hood & window integrity
- Check motor and spindle shaft for appropriate rigidity
- Check arm and motor assy for proper return
- Test front control panel for all functions
- Check and adjust handle arm tension clamp
- Check blade motor for proper operation
- Perform a cut test and monitor machine for proper operation
- Clean cutting bay
- Check coolant tank and lines for abrasive residue
- Check and lubricate pillow block for proper blade rotation
- Clean and lubricate vise beds
- Check motor and spindle bearing for proper operation
- Check pump for proper operation
- Check coolant nozzle integrity
- Replace hood gas springs



Wear Parts Kit

The parts listed below are included with every Preventative Maintenance visit

Delta Manual Wear Parts Kit 102213SRV

Includes

Motor Drive Belt

Coolant

Hood Gas Spring

Water Nozzle

Clamp

Spark Pad

Blade Spindle

Spindle Bearings



BUEHLER Worldwide Headquarters
North America-South America Offices
41 Waukegan Road
Lake Bluff, Illinois 60044-1699 USA
P: 800 BUEHLER (800-283-4537)
P: (847) 295-6500
www.buehler.com | info@buehler.com

European Headquarters
BUEHLER Germany
info.eu@buehler.com

BUEHLER France
info.eu@buehler.com

BUEHLER United Kingdom
info.eu@buehler.com

BUEHLER China
info.cn@buehler.com

BUEHLER Japan
info.japan@buehler.com

BUEHLER Asia-Pacific
info.asia@buehler.com